

— MIRAGE 2.0 SWING MOUNT PACKAGE · MSP3500-2.0 / MSP3550-2.0-CB

Mirage 2.0

Swing Mount Operator System

Swing mount operator with the 4000-2.0 hydraulic chair (steel or cast aluminum base), 3500-2.0 rotating unit and 3550-2.0 rear radius assistant's arm.



Specifications

01 MODELS

Model	Description
MSP3500-2.0	Mirage 2.0 Swing Mount Operator System with Steel Base
MSP3550-2.0-CB	Mirage 2.0 Swing Mount Operator System with Cast Aluminum Base

02 PACKAGE INCLUDES

Part	Description
4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control
3500-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/ dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system, with touch pad control
3550-2.0	Rear radius asepsis assistant's arm system, Hve and saliva ejector, water quick disconnect and solids collector, 3 Way Syringe

Part	Description
4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminium base

03 PACKING DIMENSIONS

Chair: 44 x 29 x 33 340 lbs.

4000-2.0 CB Cushion: 51 x 24 x 10 35 lbs.

5500: 36 x 28 x 39 100 lbs.

3550: 29 x 12 x 24 25 lbs.

LP100: 45 x 26 x 44 30 lbs.

04 OPTIONS AND PARTS

Part	Description
3500-2.0	Mirage swing mount unit with chair bracket (for TPC chair only)
3550-2.0	Mirage rear swing mount suction package (for TPC chair only)
LP-100	Mirage swing mount light post (ordered w/unit)
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)
L700-LED	Vivid LED Operatory Light
L795-LED	Vivid LED Ceiling Mounted Operatory Light
L550-LED	Lustrous LED Operatory Light
L7601-LED	Luminous LED Operatory Light
P2335+TB	Deluxe autoclavable Hve Long valve with tubing
WIC800	Standard flat screen monitor bracket
WIC840	Mirage flat screen monitor bracket with 16" extension
M017	17" LCD multi media flat screen monitor
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply
53016-L	Disposable Solids Collector Trap 2 1/8" (pack of 100)
EMB-901	Built in Electric Motor System
AS-1101	Mirage Assistants stool
DR1102	Mirage Doctor's Stool

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.