

— MIRAGE 2.0 OPERATORY PACKAGE · MP2000-L550LED-2.0-CB / MP2015-L550LED-2.0-CB

Mirage

Chair Mounted Operatory System, Cast Aluminum Base

Complete operatory with the 4000-2.0-CB hydraulic chair, 2015-2.0 delivery unit and L550-LED Lustrous light. Available with cuspidor or with assistant's instrumentation.



Specifications

01 MODELS

Model	Description
MP2000-L550LED-2.0-CB	Mirage 2.0 Operatory Package with Cuspidor
MP2015-L550LED-2.0-CB	Mirage 2.0 Operatory Package with Assistant's Instrumentation

02 PACKAGE INCLUDES

Part	Description
4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminium base
2015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, Hve and saliva ejector, water / air quick disconnects, solids collector clean water system Touch pad Control on Both Dr. & Assistant instrument L550-LED Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux composite cure mode, motion sensor, manual toggle on/off
2000-C	Porcelain cuspidor with timed bowl rinse and cup filler

03 PACKING DIMENSIONS

Chair: 44 x 29 x 33 340 lbs.

Cushion: 51 x 24 x 10 35 lbs.

Unit: 39 x 29 x 23 60 lbs.

Assist PKG: 43 x 25 x 12 40 lbs.

Dental Light Arm L600: 34 x 11 x 5 30 lbs.

Cuspidor 26 x 18 x 10 20 lbs.

04 OPTIONS AND PARTS

Part	Description
TPC3WAY-PRT	Additional autoclavable 3 way syringe with tubing
P2335+TB	Deluxe autoclavable Hve Long valve with tubing
WIC800	Standard flat screen monitor bracket
WIC840	Mirage flat screen monitor bracket with 16" extension Arm
M017	17" LCD multi media flat screen monitor
A800N-1	Advance 800 piezo scaler – built in
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply
ALED55BN-A	LED 5W curing light – built in to dental unit
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100)
TPE-2011	Touch pad Control - Chair & Electric Motor System touch pad
EMB-901	Built in Electric Motor System
XV100	Super thin X-ray film viewer
XV150-R	X-ray film viewer bracket (curved) for TPC Mirage unit only
AS-1101	Mirage Assistant's stool
DR-1102	Mirage Doctor's stool

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.