

Mirage

Panel & Cabinet Mount Systems

Panel mount rear delivery (MP-2600), rear cabinet mount delivery with assistant's package (MR-4500 L/R), and wall / cabinet mount assistant's package (MP-2601).



01 FEATURES

- ✓ Panel mount, cabinet mount and wall mount options
- ✓ Vacuum canister with removable screen
- ✓ Left / right assistant's package available

Specifications

02 MODELS

Model	Description
MP-2600	Mirage Panel Mount Rear Delivery System (2 HP manual select, 3-way syringe, 3 position HP holder, utility valves)
MR-4500 L/R	Mirage Rear Cabinet Mount Delivery System with assistant's package (2 HVE, 1 SE, 3 position holder)
MP-2601	Mirage Wall & Cabinet Mount assistant's package (2 HVE, 1 SE, 3 position holder)

03 MP-2600 FEATURES

- ✓ 2 handpiece tubing positions
- ✓ 3 position HP holder
- ✓ TPC press ring style syringe, with adjustable air & water controls
- ✓ Individual drive air & water coolant controls
- ✓ Foot control
- ✓ Air / water master control with regulator
- ✓ Control dimensions 8 1/2" x 1 3/4" x 5"
- ✓ Cut-out dimensions 7 1/2" x 11 1/2"

04 MR-4500 L/R FEATURES

- ✓ Asepsis 3 handpiece automatic selection with bottle system
- ✓ 3 way syringe with quick change autoclavable tip
- ✓ Master on/off toggle
- ✓ Air / water pressure adjustment for each handpiece
- ✓ Pressure gauge
- ✓ Non-retracting water / air coolant adjustment and purge
- ✓ Asepsis handpiece tubing
- ✓ Self contained bottle system
- ✓ Stainless steel tray
- ✓ Vacuum canister with removable screen
- ✓ 4 position holder, 1 autoclavable HVE and saliva ejector with tubing
- ✓ Variable speed foot control

05 MP-2601 FEATURES

- ✓ 2 HVE suction valves
- ✓ 1 SE suction valve
- ✓ 3 position suction holder
- ✓ Vacuum canister with removable screen

06 SPECIFICATIONS

MP-2600 packing dimensions	9" x 17" x 15" / 9 lbs
MP-2601 packing dimensions	8" x 17" x 15" / 5 lbs
MR-4500 L/R packing dimensions	19" x 15" x 20" / 30 lbs

07 OPTIONS AND PARTS

Part	Description
H7660	Illumination system, ISO C (6 pin): tubing, control box
H7000	Illumination system, ISO B (5 hole): tubing, control box, bulb
A800N	Advance 800 piezo scaler, built in

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.