

— ORTHODONTIC MOBILE DELIVERY CABINET · OMC2375

OMC2375 Orthodontic Mobile Delivery Cabinet

Orthodontic mobile delivery cabinet for centrally plumbed operatories, with a pullout writing surface, instrument drawers, storage drawers and oversized locking casters.



01 FEATURES

- ✓ Pullout writing surface
- ✓ Two instrument drawers, one medium and two large storage drawers
- ✓ Oversized locking dual-caster wheels
- ✓ Optional hard surface stone top (Artic, White, Cream, Bamboo)
- ✓ Optional OMC 2200 monitor bracket

Specifications

02 MODELS

Model	Description
OMC2375	Orthodontic Mobile Delivery Cabinet

03 STANDARD CABINET FEATURES

- ✓ Five easy to remove and clean full extension drawers and one pull-out writing surface
- ✓ High quality hardware throughout
- ✓ Oversized caster wheels with locks
- ✓ Four stationary handles for easy cart positioning
- ✓ A removable counter-top plastic tray
- ✓ Standard 12" x 12" junction box
- ✓ 7' length of umbilical

04 STANDARD CONTROL FEATURES

- ✓ 3 HP automatic selection
- ✓ 3-way syringe
- ✓ Vacuum package with saliva ejector and HVE
- ✓ Self contained water bottle system
- ✓ Variable speed disc-type foot control
- ✓ Water coolant on/off switch with flow control
- ✓ 1 year warranty
- ✓ Includes utility center + junction box
- ✓ Requires central air and vacuum connections

05 SPECIFICATIONS

Packing dimensions 41" x 28" x 32" / 140 lbs

Rear compartment 4" W x 12" D x 20" H

06 OPTIONS AND PARTS

Part	Description
A800N-1	Advance 800N Piezo Scaler - Built in
ALED55BN-A	LED 55BN 5W LED Curing Light -Built in
2375-UMB-FT	Umbilical Per Foot
OMC2200	LCD Monitor Bracket W/ Post
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic
OMC-ST82	Optional Stone Surface on OMC 2375, Color White
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.

SOLID STONE SURFACE COLORS



Artic
OMC-ST81



White
OMC-ST82



Cream
OMC-ST85



Bamboo
OMC-ST86