

— ORTHODONTIC MOBILE DELIVERY CABINET

OMC 2375CV

A self-contained orthodontic delivery cabinet with five full-extension drawers, a pull-out writing surface and on-board air, water and suction. Oversized dual casters roll easily on every floor, including carpet.

OMC2375CV-SL shown



5 drawers

Full extension, easy to remove and clean

3 HP

Automatic handpiece selection

4 stone tops

Optional hard-surface counters

1 year


Limited warranty



Rear compartment. Houses the dual compressor, 4" W × 12" D × 20" H.



Pull-out writing surface. Shown with optional solid surface.



Optional monitor bracket. LCD bracket with post, [OMC2200](#)

Specifications

TPC

01 STANDARD CABINET FEATURES

- ✓ Five full-extension drawers and one pull-out writing surface, all easy to remove and clean
- ✓ High-quality hardware throughout
- ✓ Oversized caster wheels with locks
- ✓ Four stationary handles for easy positioning
- ✓ Removable plastic countertop tray
- ✓ Large open space with front access for accessories

03 MODELS

Included	CV-SL	CV
All standard features	✓	✓
Built-in piezo scaler	✓	—
Built-in LED curing light	✓	—
1-year warranty	✓	✓

02 STANDARD CONTROL FEATURES

- ✓ 3 HP with automatic selection
- ✓ 3-way syringe
- ✓ Vacuum package with saliva ejector and HVE
- ✓ Self-contained water bottle system
- ✓ Variable-speed disc-type foot control
- ✓ Water coolant on/off with flow control
- ✓ Self-contained dual compressor

04 DIMENSIONS

Top	26.5" W × 22" D, angled corner
Rear compartment	4" W × 12" D × 20" H


05 PACKING

Carton	L × W × H in	Weight
OMC2375CV-SL	32 × 28 × 41	180 lb

06 OPTIONS

OMC2200 LCD monitor bracket with post

OMC-ST81 to ST86 Optional hard-surface stone countertop, in four colors



Artic
OMC-ST81



White
OMC-ST82



Cream
OMC-ST85



Bamboo
OMC-ST86

Options are not factory-installed unless specified on the purchase order. Stone colors are approximate.

Self-contained, so it goes where you need it. On-board dual compressor, water bottle system and suction mean no plumbing hookups: just roll it chairside and plug it in.