

— ORTHO 1200 ORTHODONTIC PACKAGE · ORP1200-550LED-2375 / CV / CV-SL

Ortho 1200 Fixed Base Orthodontic Package

Orthodontic package with the OR1200 fixed base chair, OMC2375 mobile delivery cabinet and L550-LED Lustrous light, available centrally plumbed or self-contained.



01 FEATURES

- ✓ OR1200 fixed base chair with last position memory
- ✓ OMC2375 mobile delivery cabinet
- ✓ L550-LED Lustrous operatory light
- ✓ Central plumbing or self contained (compressor + vacuum)
- ✓ 1 year warranty

Specifications

02 MODELS

Model	Description
ORP1200-550LED-2375	Ortho 1200 Orthodontic Package, OMC2375 cabinet
ORP1200-550LED-2375CV	Ortho 1200 Orthodontic Package, OMC2375CV: cabinet + compressor + vacuum (self contained)
ORP1200-550LED-2375CV-SL	Ortho 1200 Orthodontic Package, OMC2375CV-SL: cabinet + compressor + vacuum + LED curing light + scaler (self contained)

03 PACKAGE INCLUDES

Part	Description
OR1200	Fixed base electromechanical patient chair with wide back and full lumbar support, double articulating headrest with sliding height adjustment, one armrest included, foot control and last position memory
OMC-2375	Orthodontic mobile delivery cabinet with a pullout writing surface, two instrument drawers, one medium and two large storage drawers that fully extend and quickly remove for easy cleaning, and convenient access from the front. Oversized locking dual-casters.
L550-LED	Lustrous LED post mounted light with dual intensity, 8,000 – 30,000 Lux, dual axis light head rotation

04 CHAIR FEATURES

- ✓ Independent backrest
- ✓ Articulating dual axis headrest
- ✓ Fixed base
- ✓ Built in fuse / surge protection
- ✓ Last position memory foot control
- ✓ Short pump cover
- ✓ Electric motor on backrest

05 CABINET FEATURES

- ✓ Five easy to remove and clean full extension drawers and one pull-out writing surface
- ✓ High quality hardware throughout
- ✓ Oversized caster wheels with locks
- ✓ A removable counter-top plastic tray
- ✓ Standard 12" x 12" junction box
- ✓ 7' length of umbilical

06 LIGHT FEATURES

- ✓ Composite cure mode
- ✓ 24V DC power supply
- ✓ Bulb type: 4 white LEDs, 2 yellow
- ✓ Light intensity (Lux): high 30,000, low 8,000
- ✓ Focal length: 27" (70 cm)
- ✓ Illumination area: 100mm x 180mm
- ✓ Sensor mode, manual toggle on/off

07 CONTROL FEATURES

- ✓ 3 HP automatic selection
- ✓ 3-way syringe
- ✓ Vacuum package with saliva ejector and HVE
- ✓ Self contained water bottle system
- ✓ Variable speed disc-type foot control
- ✓ Water coolant on/off switch with flow control
- ✓ 1 year warranty
- ✓ Requires central air / vacuum systems. Includes utility center and junction box.
- ✓ Optional: self contained compressor and vacuum (no plumbing required)

08 PACKING DIMENSIONS

Chair: 44 x 29 x 33 / 340 lbs.

Cushion: 51 x 24 x 10 / 35 lbs.

Light arm L600: 26 x 14 x 8 / 30 lbs.; light head L550: 17 x 15 x 7

Cabinet: 32 x 28 x 41 / 140 lbs.

09 OPTIONS AND PARTS

Part	Description
WIC840	Mirage flat screen monitor bracket with 16" extension
M017	17" LCD multi media flat screen monitor
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply
OMC2200	LCD monitor bracket
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic
OMC-ST82	Optional Stone Surface on OMC 2375, Color White
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo
L2001	Laguna L2000-2.0 Chair Bracket
AS-5047	Laguna Assistants stool
DR-5045	Laguna Doctor's stool
PTH-101	Pony Tail Headrest

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.