

Ortho 1200 Fixed Base Orthodontic Patient Chair

A fixed base orthodontic patient chair with an independent backrest, articulating dual axis headrest, steel plate base, foot control and last position memory.



01 FEATURES

- ✓ Independent backrest
- ✓ Articulating dual axis headrest
- ✓ Steel plate base
- ✓ Fixed base chair
- ✓ Foot control and last position memory

Specifications

02 SPECIFICATIONS

Model	OR1200
Base	Steel plate, fixed base height
Backrest	Independent backrest
Headrest	Double articulating, optional pony tail headrest
Controls	Electronic foot control; one preset exit position
Chair size (L x W x H, in)	44 x 29 x 33

Model OR1200

03 PRODUCT INCLUDES

- ✓ Fixed base height
- ✓ One preset exit position
- ✓ Electronic foot control
- ✓ Independent backrest
- ✓ Toe guard protection to prevent scuffing
- ✓ Double articulating headrest, steel plate base
- ✓ Seamless upholstery

04 OPTIONS

PC-5000-DU	Upcharge: dental chair upholstery, wide back, LK vinyl (specify color and product code)
PC-5000-Vinyl-UC	Upcharge: non-standard dental chair upholstery, wide back, vinyl
PC-5000-ULTR-UC	Upcharge: Ultra leather dental chair upholstery, wide back
L2001	OR-1200 orthodontic chair bracket
PTH 101	Pony tail headrest
M42408C	Right or left side armrest with bracket

Options are not factory-installed unless specified on the purchase order. Additional charges apply. Contact TPC for pricing.