

— AIR POLISHING SYSTEM · P555 / P557

Air Polisher

Prophy power that simply connects to handpiece tubing: no installation cost, no extra foot control, and no external air or water lines.

P557 · P555



360°

Swivel nozzles, two included

134°C

Autoclavable nozzle

130g

Prophy powder included

1 year

Limited warranty



Connects in seconds. Runs from standard 4-hole handpiece tubing.



Color-coded covers. Three translucent covers included with the P555.

Specifications



01 PACKAGE INCLUDES

Included	P555-4H	P557-4H
130 g prophy powder	✓	✓
3 color-coded powder holder covers	✓	—
2 nozzles, 360° swivel	✓	✓
Cleaning tool	✓	✓
O-ring set	✓	✓
1-year warranty	✓	✓

02 DESIGN FEATURES

- ✓ **EZ nozzle.** Length, angle and curve reach anterior and posterior teeth
 - ✓ **Quick connection.** Pull back the release ring to remove the nozzle
 - ✓ **Cleaning power.** 130 g jet polishing powder
- ✓ **Ergo comfort.** Notched grip fits large and small hands
 - ✓ **No installation cost.** No extra foot control or utility lines
 - ✓ **Lightweight and flexible** for easy access

03 OPTIONS

P560	130 g prophy powder
P660	Nozzle for P555
P665	Nozzle for P557

04 PACKING DIMENSIONS

Carton	L × W × H in	Weight
Single unit	9 × 8 × 3	2 lb
Case of 10	13 × 13 × 10	12 lb

Fully autoclavable nozzle. The handpiece nozzle withstands sterilization up to 134 °C / 273 °F for reliable infection control between patients.