

— MOBILE VACUUM SYSTEM

PC 2530

A portable, self-contained vacuum for temporary locations, or wherever installing a central vacuum isn't feasible. It also makes a dependable backup.

PC2530 shown

**75**cfm

Airflow

7/8HP

Variable-speed motor

70dB

Sound level

1year

Limited warranty

STANDARD FEATURES

- ✓ Moves with ease on locking casters
- ✓ Simple maintenance, 2 L canister
- ✓ Compact metal case, 14" × 12" footprint
- ✓ Variable speed control
- ✓ Perfect backup for central vacuum
- ✓ 115 VAC or 230 VAC models

Specifications



01 PACKAGE INCLUDES

- ✓ 2 HVE valves and 1 SE valve
- ✓ 2 SE suction tips
- ✓ Solids trap
- ✓ Drain line
- ✓ Power cord
- ✓ 1-year warranty

02 VACUUM SPECIFICATIONS

Motor	7/8 HP, variable speed
Electrical	115 VAC or 230 VAC, 50/60 Hz
Current	4.0 A
Vacuum	80" of water (6" Hg)
Airflow	75 cfm
Canister	2 liter
Sound	70 dB
Construction	Metal case
Dimensions	14" L × 12" W × 31" H
Shipping	Freight

03 PACKING DIMENSIONS

Carton	L × W × H in	Weight
PC2530	17 × 25 × 37	95 lb

04 IDEAL FOR

- ✓ Temporary or satellite locations
- ✓ Operatories without central vacuum
- ✓ Backup during service or downtime
- ✓ Pairing with portable delivery systems

Pairs with TPC portable delivery. Add a PC2530 to a PC2635 or PC2630 system for full HVE suction anywhere you set up.