

— PORTABLE DELIVERY SYSTEM · PC2630 / PC2930

PC 2630

A fully self-contained air dental system that packs into a rugged transport case. Quick to set up and operate, and ready for travel. PC2930 adds a built-in piezo scaler and LED curing light.

PC2630 shown



580_w

Oil-less compressor, 7 L

2 models

PC2930 adds scaler and curing light

60 lb

Packed weight (PC2930: 65 lb)

1 year

Limited warranty



Rolling transport case. Pull-out handle and wheels for travel.



Built-in scaler and curing light. Standard on [PC2930](#)

Specifications



01 PACKAGE INCLUDES

Included	PC2630	PC2930
4-hole or 2-hole instrument tubing (specify)	✓	✓
3-way syringe	✓	✓
Saliva ejector	✓	✓
Self-contained water bottle system	✓	✓
Self-contained oil-less compressor (580W / 7L)	✓	✓
Foot control	✓	✓
Air venturi vacuum system	✓	✓
Transport case with pull-out handle	✓	✓
Built-in piezo scaler	–	✓
Built-in LED curing light	–	✓
1-year warranty	✓	✓

02 SYSTEM SPECIFICATIONS

Type	Self-contained air system
Compressor	Oil-less, 580W / 7L
Vacuum	Air venturi
Water	Self-contained bottle
Case	Transport case, pull-out handle
Handpieces	Tornado available on request
Warranty	1-year limited

03 OPTIONS

A800N-1	Advance 800 piezo scaler, built in
INSTALL-1	Labor to install built-in piezo scaler with additional holder on PC2630
53012-INSTL	Female quick disconnect, factory installed

04 PACKING DIMENSIONS

Carton	L × W × H in	Weight
PC2630	16 × 20 × 28	60 lb
PC2930	16 × 20 × 28	65 lb

Options are not factory-installed unless specified on the purchase order.

Ready for travel. Everything packs into one rugged case with a pull-out handle and wheels. Tornado handpieces are available on request; ask for details.