

— PORTABLE DELIVERY SYSTEM

PC 2635

A lightweight, self-contained delivery system with its own air and water supply. Ready for in-office use, nursing homes, hospitals or anywhere off-site dentistry is needed.

PC2635 shown

**580_w**

Oil-less compressor, 7 L

4 / 2_{hole}Instrument tubing
(specify)**75_{lb}**

Packed shipping weight

1_{year}

Limited warranty

STANDARD FEATURES

- ✓ Self-contained oil-less compressor
- ✓ Self-contained water bottle system
- ✓ 3-way syringe and foot control
- ✓ Automatic handpiece flush system
- ✓ Built-in carrying handle
- ✓ Tornado high- and slow-speed handpieces available

Specifications

TPC

01 PACKAGE INCLUDES

- ✓ 4-hole or 2-hole instrument tubing (specify)
- ✓ Self-contained oil-less compressor (580W / 7L)
- ✓ Automatic handpiece flush system
- ✓ Self-contained water bottle system
- ✓ 3-way syringe
- ✓ Foot control
- ✓ Built-in carrying handle
- ✓ 1-year warranty

02 SYSTEM SPECIFICATIONS

Type	Portable, self-contained
Compressor	Oil-less, 580W / 7L
Air and water	Supplied on board
Instrument tubing	4-hole or 2-hole
Handpiece care	Automatic flush system
Handpieces	Tornado high and slow speed, optional
Carrying	Built-in handle
Warranty	1-year limited

03 PACKING DIMENSIONS

Carton	L × W × H in	Weight
PC2635	23 × 17 × 22	75 lb

04 BUILT FOR CARE ON THE GO

In-office

A second chairside unit when you need one.

Nursing homes

Bring treatment to residents.

Hospitals

Bedside procedures without plumbing.

Off-site

Outreach, missions and mobile clinics.

Air and water on board.

A self-contained compressor and water bottle system mean there is nothing to plumb. Pick it up by the built-in handle, plug it in and start treating.

