

Scaler Inserts

Universal and slim ultrasonic scaler inserts in 25K and 30K, with metal or plastic handles. Interchangeable with other brands of scalers.

S1025 / S1030



S2025 / S2030



S3025 / S3030



S4025 / S4030



8 models

Universal and slim, 25K and 30K

2 water

External or internal delivery

25K / 30K

Matches your scaler frequency

90 day

Warranty against defects



Black grip. Non-slip ergonomic plastic handle. [S3025](#) / [S3030](#)



Blue grip. Slim tip with the same ergonomic handle. [S4025](#) / [S4030](#)



Metal handle. External water, autoclavable stainless. [S1025](#) / [S1030](#)

Specifications

01 MODELS

Model	Freq.	Style	Water	Handle
S1025	25K	Universal	External	Metal
S1030	30K	Universal	External	Metal
S2025	25K	Universal slim	External	Metal
S2030	30K	Universal slim	External	Metal
S3025	25K	Universal	Internal	Black plastic
S3030	30K	Universal	Internal	Black plastic
S4025	25K	Universal slim	Internal	Blue plastic
S4030	30K	Universal slim	Internal	Blue plastic

02 INSERT STYLES



S1025 / S1030 Universal, metal



S2025 / S2030 Slim, metal



S1025 / S1030 Universal, metal



S2025 / S2030 Slim, metal



S3030 Universal, black plastic



S3025 Universal, black plastic



S4025 Slim, blue plastic



S4030 Slim, blue plastic

03 FEATURES

- Interchangeable with other brands of scalers
- External water (metal handle) or internal water delivery (plastic handle)
- Slim tip designs available
- 90 day warranty against defects

04 PLASTIC HANDLE

- Water delivered directly to the insert tip
- Large, round, non-slip grip for tactile sensitivity and less hand fatigue
- Better water management for greater patient comfort

Fits the scaler you already own.

25K and 30K inserts interchange with other brands, so you can stock a single insert line for every operator.