

LED X-RAY VIEWER · XV100

Super Thin LED X-Ray Viewer (XV100)

0.5" super thin LED viewer for countertop, wall or unit head mounting, with instant illumination and no grey in the viewing area.



01 FEATURES

- ✓ Viewing area 8.5" H x 11" W x 0.5" thick
- ✓ Countertop, wall or unit head mounted
- ✓ AC power source with sealed LED lights
- ✓ Instant illumination, no grey in viewing area
- ✓ Accepts panoramic, cephalometric and regular size film
- ✓ 1 year warranty

Specifications

02 MODELS

Model	Description
XV100	Super Thin X-ray Viewer

03 OPTIONS AND PARTS

Part	Description
XV150-R	X-ray viewer bracket (round cut-out), mounts to Mirage unit head (round arm)
XV150-F	X-ray viewer bracket (flat cut-out), mounts to unit head (flat arm)
XV130	Replacement LED light strip
XV100-PSW	Viewer replacement power supply
XV100-CLIP	Viewer replacement clip (3 pack)

04 PACKING DIMENSIONS

17" x 11" x 2" / 3 lbs.

(10/cs.) 17" x 11" x 19" / 28 lbs.