

— TMC MOBILE CABINETS · TMC-180 / TMC-140 / TMC-160 / TMC-120

TMC Doctors & Assistants Mobile Cabinets

Mobile cabinets with high pressure laminate exteriors and ball bearing glides: utility, doctor's and two assistant's models, in white or light grey.



01 FEATURES

- ✓ Utility, doctor's and assistant's models
- ✓ High pressure laminate exterior, white interior
- ✓ Ball bearing glides and durable casters
- ✓ White or light grey
- ✓ 1 year warranty

Specifications

02 MODELS

Description	Light grey	White
Utility Mobile Cabinet	TMC-180-G	TMC-180-W
Doctor's Mobile Cabinet	TMC-140-G	TMC-140-W
Assistant's "Alabama" Mobile Cabinet	TMC-160-G	TMC-160-W
Assistant's "North Carolina" Mobile Cabinet	TMC-120-G	TMC-120-W
Assistant's "North Carolina" Mobile Cabinet, RT pull out	—	TMC-120R-W

03 TMC-180 UTILITY CABINET

- ✓ 14.5" W x 16" D x 32" H
- ✓ Designed for ultrasonic scalers, curing lights, electrosurge, intraoral cameras, etc.
- ✓ 1 instrument drawer, 2 adjustable shelves
- ✓ Shipping dimensions: 20" x 22" x 34" / 65 lbs

04 TMC-140 DOCTOR'S CABINET

- ✓ 21.5" W x 19" D x 32" H
- ✓ 4 instrument drawers, 3 storage drawers
- ✓ Shipping dimensions: 33" x 22" x 24" / 110 lbs

05 TMC-160 ASSISTANT'S "ALABAMA" CABINET

- ✓ 21.5" W x 19" D x 32" H
- ✓ Sliding counter top with amalgamator well; side to side sliding counter
- ✓ 1 instrument drawer, 2 storage drawers
- ✓ Shipping dimensions: 33" x 22" x 24" / 110 lbs

06 TMC-120 ASSISTANT'S "NORTH CAROLINA" CABINET

- ✓ 21.5" W x 19" D x 32" H
- ✓ Solid top with pull handle; amalgamator well drop leaf door
- ✓ 3 instrument drawers, 1 storage drawer
- ✓ Shipping dimensions: 33" x 22" x 24" / 110 lbs

07 ALL CABINETS

- ✓ External surfaces of high pressure laminate with complimentary white interior surface for easy cleaning
- ✓ Smooth flowing ball bearing glides and durable casters
- ✓ Color choice: white or light grey (white is shipped if color is not indicated on the purchase order)
- ✓ Mobile cabinets ship on pallet via freight carrier
- ✓ 1 year warranty