

DENTSONIC ULTRASONIC CLEANERS · UC-400 / 1000 / 1500

UC-400 / 1000 / 1500

Powerful 40 kHz cavitation for thorough, quieter cleaning. Three large tank sizes from 3.8 to 15.8 quarts, each with a stainless steel tank and cover, digital timer, basket and lid.



UC-400
3.8 qt.



UC-1000
10.6 qt. · drain kit



UC-1500
15.8 qt. · drain kit



Drain kit
Ball valve drain on UC-1000
and UC-1500

40 kHz

Cavitation for thorough cleaning

3 sizes

3.8, 10.6 and 15.8 qt.

1-99 min

Digital LCD timer

1 year

Limited warranty

Stainless steel tank and cover. Welded drain on the UC-1000 and UC-1500, vinyl clad steel housing.

Sweep technology. Improves cleaning efficiency across the tank.

Basket and lid included. Replacement baskets: **UC400-BASKET UC1000-BASKET UC1500-BASKET**

Specifications

01 MODELS

	UC-400	UC-1000	UC-1500
Input	115V / 60Hz; 230V / 50Hz	115V / 60Hz; 230V / 50Hz	115V / 60Hz; 230V / 50Hz
Frequency	40 kHz	40 kHz	40 kHz
Overall dims (LxWxH)	10.6" x 6.3" x 9.5"	14.2" x 10.6" x 12.2"	15" x 13" x 12.5"
Tank dims (LxWxH)	9.5" x 5.5" x 4"	11.8" x 9.5" x 5.9"	13" x 11.8" x 5.9"
Tank capacity	3.8 qt. / 3.59 L	10.6 qt. / 10 L	15.8 qt. / 15 L
Timer	LCD, 1-99 min	LCD, 1-99 min	LCD, 1-99 min
Stainless steel tank	Yes	Yes	Yes
Ship weight	10 lbs.	19 lbs.	24 lbs.

02 PACKAGE INCLUDES

Package includes	UC-400	UC-1000	UC-1500
Ultrasonic cleaner	✓	✓	✓
Stainless steel basket	✓	✓	✓
Lid	✓	✓	✓
Drain kit		✓	✓
Digital timer	✓	✓	✓

03 FEATURES

- Powerful 40 kHz cavitation, quieter cleaning
- Stainless steel cover and tank with welded drain
- Durable vinyl clad steel housing
- Digital timer, 1-99 minutes
- Sweep technology for cleaning efficiency
- Basket and lid included, 1 year warranty

04 OPTIONS

UC400-BASKET	Replacement basket (3.8 qt.)
UC1000-BASKET	Replacement basket (10.6 qt.)
UC1500-BASKET	Replacement basket (15.8 qt.)