

— DENTSONIC ULTRASONIC CLEANERS · UC-250 / 300 / 600 / 750

UC Series

Powerful 35 kHz cavitation for thorough, quieter cleaning. Four tank sizes from 2.5 to 7.5 quarts, each with a stainless steel tank, digital timer, basket and lid.



UC-250
2.5 qt.



UC-300
3.0 qt.



UC-600
6.0 qt.



UC-750
7.5 qt.

35 kHz

Cavitation for thorough cleaning

4 sizes

2.5, 3.0, 6.0 and 7.5 qt.

1-30 min

Digital LCD timer

1 year

Limited warranty



UC-300 basket. Fits the 3.0 qt. tank. **UC-300**



Stainless basket. Included with the 6.0 qt. tank. **UC-600**



Glass beaker set. 500 ML x 2 and holder. **UC-600-GBS**

Specifications

01 MODELS

	UC-250	UC-300	UC-600	UC-750
Input	115V / 230V	115V / 230V	115V / 230V	115V / 230V
Frequency	50/60 Hz	50/60 Hz	50/60 Hz	50/60 Hz
Overall dims (LxWxH)	12.6" x 9.6" x 8.0"	17.5" x 10.8" x 8.1"	21.0" x 12.0" x 9.6"	21.3" x 14" x 10.9"
Tank dims (LxWxH)	9.8" x 6.0" x 3.0"	9.8" x 7.0" x 3.0"	13" x 8.0" x 4.0"	13" x 9.4" x 4.1"
Tank capacity	2.6 qt. / 2.46 L	3.3 qt. / 3.028 L	6.4 qt. / 6.056 L	8 qt. / 7.5 L
Timer	LCD, 1-30 min	LCD, 1-30 min	LCD, 1-30 min	LCD, 1-30 min
302 stainless tank	Yes	Yes	Yes	Yes
Ship weight	8 lbs.	13 lbs.	20 lbs.	27 lbs.
Shipping carton	14" x 11" x 10"	20" x 14" x 12"	23" x 16" x 14"	24" x 16" x 14"

02 PACKAGE INCLUDES

Package includes	UC-250	UC-300	UC-600	UC-750
Ultrasonic cleaner	✓	✓	✓	✓
Basket	Plastic	Stainless	Stainless	Stainless
Plastic lid	✓	✓	✓	✓
Drain kit		✓	✓	✓
Digital timer	✓	✓	✓	✓
Heater	✓	✓	✓	✓
Degas feature	✓	✓	✓	✓

03 FEATURES

- Powerful 35 kHz cavitation
- Stainless steel tank
- Scratch-resistant, easy-clean plastic housing
- Digital timer
- Sweep technology for exceptional cleaning
- Basket and lid included, 1 year warranty

04 OPTIONS

UC-600-GBS	Glass beaker set, 500 ML x 2 and holder
UC600-BSKT-P	Replacement plastic basket
UC600-BSKT-SS	Replacement stainless basket
UC750-BSKT-SS	Replacement stainless basket