

— WET-RING VACUUM · WV3 / WV5

Superb Vac

A quiet single-pump wet-ring vacuum for one to five users, built in aluminum-bronze with a stainless steel impeller. Recycler models cut water use by up to 80%.

WV3 / WV5 shown



1-5 users

Simultaneous users

1,500 L/min

Max suction (WV3: 900 L/min)

58 dB(A)

Noise level

5 year

Full product limited warranty

STANDARD FEATURES

- ✓ Durable aluminum-bronze construction
- ✓ Stainless steel impeller
- ✓ Computer-balanced, sealed motor bearings
- ✓ Easy-clean filter with clear bowl
- ✓ Water recycling saves up to 80% (R models)
- ✓ UL / cUL classified

Specifications

TPC

01 CHOOSE YOUR MODEL

Model	WV3	WV3R recycler	WV5	WV5R recycler
Simultaneous users	3	3	5	5
Suction rate	900 L/min	900 L/min	1,500 L/min	1,500 L/min
Horsepower	1.2 HP	1.2 HP	2.0 HP	2.0 HP
Voltage	110V or 220V	110V or 220V	220V	220V
Current	8 A	8 A	13 A	13 A
Water use	4 L/min	0.8 L/min	4 L/min	0.8 L/min
Recycler & air/water separator	—	Included	—	Included
Weight	62 lb	62 lb	66 lb	66 lb
Dimension (in)	17×18×20	21×23×23	17×18×20	21×23×23
Product size L×W×H	17 × 14 × 18 in, all models			
Vacuum level	8–15 inHg, factory set to 10 inHg			
Noise level	58 dB(A)			
Frequency · breaker	50/60 Hz · 20 A circuit breaker			
Connections	Inlet 3/4" · drain 1.25" · city water 3/8" · air exhaust 2 3/8"			
Inlet water pressure	10 psi / 0.7 bar			
Ambient temperature	2 °C – 40 °C			
Warranty	5-year full product limited warranty			

02 ACCESSORIES



WV-AWS
Single-pump
air/water separator



CP-3
Vac / Air / Water



CP-4
Vac1 / Vac2 / Air /
Water



CP-2AV
Air / Vac



CP-3AV
Air / Vac1 / Vac2



CV22
Stacking rack

03 COMPLETE THE UTILITY ROOM



**Superb Air
DC6 Series**
Oil-less compressor,
1–9 users



Superb Vac Twin
WV6 / WV10
For 6–10 users



Stack them
CV22
Compressor over
vacuum on one rack