

— TWIN WET-RING VACUUM · WV6 / WV10

Superb Vac Twin

Quiet twin-pump suction for small to medium practices. Run one or both motors to match demand, with a second pump for backup that minimizes downtime and costly interruptions.

WV6(R) / WV10(R) shown



1-10 users

One or two motors run as needed

3,000 L/min

Max suction (WV6: 1,800 L/min)

65 dB(A)

Noise level

5 year

Full product limited warranty

STANDARD FEATURES

- ✓ Twin pumps with backup protection
- ✓ Adjustable vacuum, 8-15 inHg
- ✓ Aluminum-bronze, stainless steel impeller
- ✓ Water recycling saves up to 80% (R models)
- ✓ Air/water separator included
- ✓ Low-voltage control panel connections

Specifications

TPC

01 CHOOSE YOUR MODEL

Model	WV6	WV6R recycler	WV10	WV10R recycler
Simultaneous users	6	6	10	10
Suction rate	1,800 L/min	1,800 L/min	3,000 L/min	3,000 L/min
Horsepower	2.4 HP (2×1.2)	2.4 HP (2×1.2)	4.0 HP (2×2)	4.0 HP (2×2)
Current	8 A × 2	8 A × 2	13 A × 2	13 A × 2
Water use	6 L/min	1.01 L/min	6 L/min	1.01 L/min
Water recycler	—	Included	—	Included
Product weight	176 lb	176 lb	198 lb	198 lb
Shipping weight	188 lb	188 lb	208 lb	208 lb
Product size L×W×H	28 × 17 × 18 in, all models			
Shipping dimensions	32 × 25 × 24 in, all models			
Vacuum level	Adjustable 8–15 inHg, factory set to 10 inHg			
Noise level	65 dB(A)			
Power	220V · 50/60 Hz · 20 A × 2 circuit breakers			
Connections	Inlet 1" · drain 1.25" · city water 3/8" · air exhaust 2 3/8"			
Inlet water pressure	10 psi / 0.7 bar			
Ambient temperature	2 °C – 40 °C			
Warranty	5-year full product limited warranty			

02 ACCESSORIES



CP-3AV
3-button control
panel: Vac1 / Vac2 /
Air



CP-4
4-button control
panel: Vac1 / Vac2 /
Air / Water



CV22
Stacking rack, 33½ ×
22 × 26"



VC-3001 / VC-3002
Disposable vacuum canisters
Small and large, with mesh screen

03 COMPLETE THE UTILITY ROOM



Superb Air
DC6 Series
Oil-less compressor,
1–9 users



Superb Vac
WV3 / WV5
Single pump, 1–5
users



Stack them
CV22
Compressor over
vacuum on one rack