



Features

- Quiet fan operation
- Multi-arch whitening application
- Seven high intensity LED's
- Digital display with audio feedback
- High power Blue and UV light emitting diodes
- Compact mobile stand for easy storage
- Adjustable bleaching time with microprocessor-controlled timer
- Smart auto-select power: 110-240 Volts AC, 50/60 Hz
- LED bleaching system includes illuminator, mobile base, adjustable arm, power cord & user manual
- 1 year warranty
- Packing dimensions: 21" x 21" x 6" / 14 lbs. gross, 11.5 lbs. net

Specifications

Broad Spectrum	430-490 nm & 395-415 nm
Light Intensity	1,500 mW/cm ²
Voltage Output	110-240 V, 50/60 Hz
Dimensions	7.25" (L) x 8" (W) x 5.5" (H)
Weight	10 lbs.