



Features:

- Light may be adjusted diagonally, vertically, and horizontally
- Sharply defined rectangular light pattern
- 2-Position intensity switch
- Dichroic coated reflector filters out infrared rays providing cool light
- Easy open safety shield allows aseptic cleaning
- 1 year warranty
- Packing dimensions: 32" x 16" x 13" / 25 lbs.

SPECIFICATIONS

Input Voltage:	110V, 60Hz
Output Voltage:	12V High, 10.5V Low
Power Output:	55W High, 50W Low
Bulb Type:	Tungsten Halogen Type 12V-55W
Light Output:	24,000 Lux High 12,000 Lux Low
Color Temperature:	4500° Kelvin
Focal Distance:	27"
Light Pattern:	3" x 7" @ 27"