



15" Deluxe LCD Monitor

with Built-in X-ray Viewer

Active-Matrix TFT LCD Monitor

- 1024 x 768 resolution
- 500:1 high contrast ratio
- 160-degree viewing angle

Mirror Image Function

- Enables viewing from opposite perspective
- Better communication between Doctor and patient

Built-In Stereo Sound System

- Dual built-in stereo speakers



Anti-Flash and Anti-Glisten Intensified Glass Sheet

- Protects thin film transistor LCD panel
- Longer LCD life

Signal Hub Device

- Monitor for: pc, satellite tv, cable tv, dvd, digital camera, video camera, intraoral camera

Built-In Buffer Storage for Images

- Built-in storage for immediate comparison before and after treatment



Built-In Xray Viewer

- Dry or wet film
- Enhanced 50:1 enlargement
- Images easily stored into computer



Versatile One-Touch Remote Control

- Built-in infrared remote mouse, decimal channel selector, laser pointer, mirror imaging, function selector, adjustment for contrast, luminance and color

Wall, Post or Countertop Mounted

- Multiple monitors may be networked
- VESA wall mountable



MTFT15 with optional WIC840



WIC800




MTFT15 with optional WIC825

Packing dimensions: 19" x 10" x 20" / 15 lbs.



15" and 17" LCD Monitors

Multimedia Monitor and Brackets



Standard LCD Monitor Bracket

- Easily position the angle of your monitor for maximum comfort and productivity
- 4-way pivot adjustment



Deluxe LCD Monitor Bracket

- 16" extension arm
- vertical and horizontal travel
- 4-way pivot adjustment



Mirage Monitor Bracket

- 16" extension arm
- horizontal travel
- 4-way pivot adjustment



Standard and Deluxe LCD Monitor Bracket may be wall mounted

Features:

- Cable TV, DVD, Intraoral Camera READY!!
- Full 15" viewable TFT Active Matrix flat panel display
- Auto adjustment function
- Built-in dual speakers
- On screen display controls (OSDC)
- Recommended resolution: 1024 x 768
- Universal 100-240 VAC, 50/60 Hz
- Includes: AC adapter, remote control, video cable (RCA), s-video cable, PC cable and monitor
- Connections: cable TV-in, video-in, s-video-in, VGA-in, audio-in, DC-in
- 1 year warranty
- Packing dimensions: 15" Monitor 24" x 19" x 6" / 15 lbs.
17" Monitor 24" x 19" x 6" / 20 lbs.