


Model LED60



Features

- Low battery signal – visual and audible alarm
- Easy to read – LCD display
- High light intensity – 5 to 10 second cure
- No down time
- Flexibility – corded and cordless light cure
- Quiet fan – extends life of LED
- Sleep mode – conserves power, saves battery life



Specifications and Accessories

Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz
 Dimensions: 10"L x 2"W x 8"H
 Weight: 24 ounces

Handpiece

Light source: 5W light emitting diode
 Wavelength: 430-490 nm
 Light Intensity: 1,400 mW/cm²
 Battery: Lithium ion
 Warranty: 1 year



Description

LED 60N Cordless Curing Light

Model

LED-60N

Package includes the following:

LED 60N Cordless Curing Light	✓
Charging Base	✓
Protective Eye-shield	✓
Countertop Stand with Power Indicator	✓
8mm Light Guide	✓
Built-in power transformer	✓
AC power supply	✓

Options:

P108	Replacement light guide (8mm)
P208	Replacement protective eye-shield
LED 60N-BATT	Replacement rechargeable battery for LED 60
LM200	Light Meter

Standard Accessories

- LED 60N Handpiece
- Charging base
- 8mm light guide
- Eye protector
- Lithium ion battery
- Power supply
- User's Manual

Shipping Dimensions: 10"L x 7"W x 6"H / 3 lbs.
 (10/cs.) 21"L x 14"W x 10"H / 19 lbs.