



Model #L800

Model #L900-1
Ceiling Mounted
Lucent Operatory
Light

Features:

- The contoured handles are conveniently located for operator positioning and asepsis protection
- The bulb is easily removed for replacement
- The dichroic-coated reflector filters out the infrared rays and dissipates heat allowing only cool, white light to illuminate oral cavity
- The two-way intensity switch allows for correct oral cavity illumination and a reduction in eye fatigue
- Optional: motion sensor
- 1 year warranty
- Packing dimensions: 33" x 27" x 13" / 35 lbs.

SPECIFICATIONS

| | |
|--------------------|-----------------------------------|
| Input Voltage: | 110/220V, 60/50Hz |
| Output Voltage: | 12V High, 10.5V Low |
| Power Output: | 55W High, 50W Low |
| Bulb Type: | Tungsten Halogen Type 12V-55W |
| Light Output: | 29,000 Lux High 17,000 Lux Low |
| Color Temperature: | 4500° Kelvin |
| Focal Distance: | 28" |
| Light Pattern: | 3.6" x 7.1" @ 28" |