



Mirage

Orthodontic Package

Model 3000 Orthodontic Chair

- Cast aluminum with steel reinforcement
- Hydraulic base & back
- Synchronized back & seat movement
- Adjustable backrest speed
- 2 programmable pre-positions
- Trendelenburg position
- 60° horizontal swivel base
- Double articulating headrest
- Trip pan safety shut-off
- Automatic return/Last position presets
- Foot Control



Dental Light Specifications:

Dimmer mode: From 8,000-30,000 Lux
 Input voltage: 110V/220V, 60/50 Hz
 Output voltage: 24 Volts AC / DC
 Light intensity (Lux): 8,000 - 30,000 Lux
 Cooling System: Static Heat sink (no fan)
 Illumination area: 100mm x 180mm @ 700mm



Model #MOP3000

Model OMC2375 Orthodontic Cabinet

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface
- Lockable drawers
- High quality hardware throughout
- Oversized caster wheels
- Four stationary handles for easy cart positioning
- A removable countertop plastic tray
- Four-outlet power strip
- Color: white
- 1 Year warranty

Optional Items:

- Piezo scaler
- LED curing light
- Fiber optic handpiece control

Standard Control Features:

- Automatic selection for one high speed and one low speed handpiece
- 3-way syringe
- Vacuum package with saliva ejector and HVE
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- Master control switch
- Air pressure regulator

Packing Dimensions (L x W x H)	Weight
Chair: 54" x 34" x 31"	310 lbs.
Cushion: 51" x 9" x 23"	24 lbs.
Light: 37" x 17" x 15"	28 lbs.
Cabinet: 26" x 22" x 34"	115 lbs.