



Mirage

Swing Mount Operator System

Standard Features

- Asepsis automatic control for three handpieces
- Individual handpiece pressure adjustments
- Non-retracting water coolant with on/off
- Individual water coolant flow controls
- Air coolant flow control
- Handpiece tubing flush system
- Ultra soft handpiece tubing
- 4-hole straight gray asepsis tubing with Midwest connector
- Stainless steel tray with non-skid pad
- Disc foot control wet/dry toggle with 7' of smooth tubing
- Equipment utility center with asepsis housing
- Self contained water system w/city selector switch
- Pre-wiring for fiber optic light source
- Air-actuated counterbalanced arm lock

Assistant's Instrumentation

- Rear radius swing arm with HVE and SE
- Assistant's arm has 4-position holder for instruments
- QD water connection with flow control
- Deluxe contoured junction box
- Vacuum canister with solids collector screen

Patient Chair

- Hydraulic patient chair with wide back and full lumbar support
- Rotating armrest left or right
- Double articulating headrest with sliding height adjustment
- 30 degree swivel from left and right of center
- 3 programmable positions with automatic return and foot control



Packing Dimensions (L x W x H)	Weight
Chair: 55 x 33 x 33	350 lbs.
Cushion: 50 x 8 x 22	20 lbs.
Unit: 51 x 21 x 25	65 lbs.
46 x 39 x 9	62 lbs.
46 x 39 x 9	25 lbs.

