

# TPC

# Mirage

Swing Mount Operatory System

Model #MSP3500



Model #MSP3500 with Dental Light



Model #MSP3500 with LCD



Model #4000 with LCD

# Mirage Swing Mount Operatory System

## Model 4000 Dental Chair

### Standard Features

2 operator-programmable position, auto return and last position  
 Electronic foot control switch preposition / auto return switch  
 Fuse / surge protection  
 Toe guard protection  
 Synchronized seat and toe board tilt  
 All metal construction  
 Thin back with lumbar support  
 Double articulating headrest  
 Safety shut-off plate  
 Pivoting armrest  
 Seamless asepsis upholstery  
 Input voltage (110 V or 220 V)- specify  
 5 year warranty

## Mirage Swing Mount Delivery System

### Standard Features

Asepsis automatic control for three handpieces with syringe  
 Individual handpiece pressure adjustments  
 Non-retracting water coolant with on/off  
 Individual water coolant flow controls  
 Air coolant flow control  
 Handpiece tubing flush system  
 Ultra soft handpiece tubing  
 4 hole straight gray asepsis tubing with Midwest connector  
 Stainless steel tray with non-skid pad  
 Disc foot control wet/dry toggle with 7' of smooth tubing  
 Equipment utility center with asepsis housing  
 Vacuum canister with solids collector screen  
 Self contained water system with city/ bottle selector  
 Pre-wiring for fiber optic light source  
 Air-actuated counterbalanced arm lock

## Mirage Rear Radius Vacuum Package

### Standard Features

Water outlet with flow control  
 Air outlet with flow control  
 Vacuum canister with solids collector screen  
 Four position holder with HVE and saliva ejector  
 One autoclavable HVE with asepsis tubing  
 One autoclavable saliva ejector with tubing

## Dental Operatory Light

**L700-1** Advance Operatory Light, 12,000 lux / 24,000 lux, 4,500° Kelvin  
**L825** Lucent Operatory Light (sensor), 29,000 lux / 17,000 lux, 4,600° Kelvin  
**L800** Lucent Operatory Light (no sensor), 29,000 lux / 17,000 lux, 4,600° Kelvin  
 (Must order #LP100 Swing mount light post with all dental lights)

|                                     | 4000 (Dental Chair) | 3500 (Swing Mount Unit) | 3550 (Rear Radius Vacuum Unit) | L700 (Standard Operatory Light) | L825 (Sensor Operatory Light) | L800 (No Sensor Light) |
|-------------------------------------|---------------------|-------------------------|--------------------------------|---------------------------------|-------------------------------|------------------------|
| Standard Features (Dental Chair)    | ••••••••••••••••    |                         |                                |                                 |                               |                        |
| Standard Features (Delivery System) |                     | ••••••••••••••••        |                                |                                 |                               |                        |
| Standard Features (Vacuum Package)  |                     |                         | ••••••••••••••••               |                                 |                               |                        |
| Dental Operatory Light              |                     |                         |                                | •                               | •                             | •                      |




## Factory Options

| Factory Options  | Item Number |
|--|-------------|
| Additional 3 way syringe                                   | #TPC3WAY+T  |
| Additional HVE with tubing                                 | #P2335+TB   |
| Power optic system, ISO B (5-hole)                         | #H7000      |
| Power optic system, ISO C (6-pin)                          | #H7660      |
| Built-in piezo scaler with transformer                     | #800N       |
| Super thin x-ray viewer- unit mounted                      | #XV100/150  |
| Standard LCD monitor bracket                               | #WIC800     |
| Deluxe LCD monitor bracket                                 | #WIC825     |
| Mirage LCD monitor bracket                                 | #WIC840     |
| 15" standard multi media LCD monitor                       | #M015       |
| Plastic utility shelf- (grey, white or beige)-post mounted | #US100      |
| Touch pad-unit head mounted                                | #TP-2000    |
| Swing mount light post                                     | #LP100      |
| Swing mount monitor post                                   | #LP105      |

## Packing Dimensions

|            |                 |          |
|------------|-----------------|----------|
| Chair:     | 54" x 34" x 31" | 310 lbs. |
| Cushion:   | 51" x 9" x 23"  | 24 lbs.  |
| Unit:      | 44" x 36" x 27" | 121 lbs. |
| Rear Unit: | 21" x 20" x 16" | 23 lbs.  |

## Upholstery Selections:



DR-1102

