

Function:

This orthodontic mobile delivery cabinet offers a pullout breadboard, two instrument drawers, one medium and two large drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide a large open space for accessories with convenient access from the front. The oversized dual-caster wheels, combined with the four-stationary handles, makes positioning of the cabinet simple and easy on all floor surfaces, even carpet.




White

Cream

Artic

Bamboo

Optional hard surface stone top



Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface
- Lockable drawers
- High quality hardware throughout
- Oversized caster wheels
- Four stationary handles for easy cart positioning
- A removable countertop plastic tray
- Four-outlet power strip
- Color: white
- 1 Year warranty

Optional Items:

- Piezo scaler
- LED curing light
- Fiber optic handpiece control

Standard Control Features:

- Automatic selection for one high speed and one low speed handpiece
- 3-way syringe
- Vacuum package with saliva ejector and HVE
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- Master control switch
- Air pressure regulator
- Requires central air and vacuum connections

Packing Dimensions

Cabinet:	26" x 22" x 34"	115 lbs.
----------	-----------------	----------

