



Power Light 100





Halogen Curing Light



Specifications

- Output light intensity: 500-600 mW/cm²
- Output wavelength range: 400-500 nm
- Preset curing times from 10 to 60 seconds or continuous cure
- LCD digital display
- Maintains constant light output and curing strength
- 8 mm light guide, protective eye shield and instruction manual included
- Powerful 75 Watt bulb
- Ultra quiet fan
- 1 year warranty (bulb excluded)
- Packing dimensions: 12" x 12" x 12" / 8 lbs.

Lighting Guides

 P101	13mm straight for large facial surface
 P102	13mm x 80° for posterior occlusal curing
 P103	8mm x 60° for general curing
 P104	3mm x 90° for layer curing of posterior restoration