

Miscellaneous

Autoclavable High Vacuum Valve – Short



Part No.	Description	Retail
#2000-HVESHORT	Autoclavable high volume ejector–short	\$24.00
#P2033	Replacement HVE lever (5 pcs/Bag)	\$10.00/Bag

- Autoclavable valve
- Finger tip lever control for on/off flow
- Swivel quick disconnect for 360° rotation
- Adapts to 1/2" ID vacuum tubing
- Disconnect at swivel for simple autoclaving
- Anodized aluminum one piece valve body construction with anodized aluminum swivel

Autoclavable High Vacuum Valve – Long



Part No.	Description	Retail
#2000-HVE	Autoclavable high volume ejector–long	\$24.00
#P2033	Replacement HVE lever (5 pcs/Bag)	\$10.00/Bag

- Autoclavable valve
- Finger tip lever control for on/off flow
- Swivel quick disconnect for 360° rotation
- Adapts to 1/2" ID vacuum tubing
- Disconnect at swivel for simple autoclaving
- Anodized aluminum one piece valve body construction with anodized aluminum swivel

Autoclavable Saliva Ejector Valve



Part No.	Description	Retail
#2000-SE	Autoclavable saliva ejector valve	\$18.00
#P2036	Replacement saliva ejector lever <u>only</u> (5 pcs/Bag)	\$10.00/Bag
#P2039	Replacement saliva ejector rubber tip <u>only</u> (5 pcs/Bag)	\$9.00/Bag

- Autoclavable valve
- Finger tip lever control for on/off flow
- Swivel quick disconnect for 360 degree rotation
- Disconnect at swivel for simple autoclaving
- Anodized aluminum one piece valve body construction with anodized aluminum swivel