

TPC 3D




User Menu of Aidite Cloud for Lab

Intraoral Scanner User Guide

Everyone with a healthy and beautiful smile

Introduction


Aidite cloud is an integrated dental platform that support the workflow between clinics and labs.
Adress : <https://intraoralscan.aiditecloud.com/> Or you can access the cloud through the software.
Settings > About > Access Aiditecloud



1.4.1 User Menu of Aidite Cloud for Lab

Dashboard

In Dashboard, you can check cloud storage, order type distribution, order status.



My Account

In this page, you can edit your account information.

Dashboard

My Account

Case

Cloud Storage

Partner

Device

My Account My Organization

Account Information

Email Address Mobile Phone Name

aiditertest_lab3@aidite.com Please input aiditertest_lab3

Email Reminders

☐ Case sending ☒ Order sending ☒ Order accepted ☒ Order rejected ☒ Order cancelled ☒ Order shipped ☒ Order completed

Password Change

Change

Time Zone (Current Time: 2026-07-22 16:49:13 UTC+8)

UTC+08:00 Beijing, Chongqing

Save

Case

You can check the case data you received and preview the details.

Dashboard

My Account

Case

Cloud Storage

Partner


Device

All Pending Accepted Under design Producing Shipped Completed Rejected Canceled

Order No./ Patient Case Types Upload start... Upload end... Search Reset


Order No.	Sender	Patient	Type	Due Date	Upload Date	Status	Details	Remarks	Operation
20260718703535	aiditertest_clinic1	A B	Restoration		2026-07-18 15:51:24	Accepted			Under design
20260717536318	aiditertest_clinic1	DemoTest	Health Report		2026-07-17 14:57:25	Producing			Under design Producing Shipped
20260717184384	aiditertest_clinic1	DemoTest	Health Report		2026-07-17 13:34:26	Producing			Under design
20260717438587	aiditertest_clinic1	DemoTest	Health Report		2026-07-17 13:32:31	Accepted			Under design
20260717757528	aiditertest_clinic1	DemoTest	Health Report		2026-07-17 13:30:24	Pending			Audit


1 Total 5 10/page



Cloud Storage


In this page, you can check the case sent to the cloud by yourself and preview the details.






Partner

On this interface, you can search for other users, send partnership requests, accept/ reject partnership invitations and manage partnerships.





Device

On this screen, you can manage your bound intraoral scanners.

