


TPC 3D



Buttons Function of Scanner

Intraoral Scanner User Guide

Everyone with a healthy and beautiful smile




□ **-Primary Control Button**

Click:Start/stop scanning Long Press:Enter Remote Control Interface

□ **-Motion-sensing Button**

Click:Automatic Switch to Next Scanning Object Long Press:Enter Motion Control Mode to rotate the model,Release the button to exit .



Back

□ **-Secondary Control Button**

Click:Start/stop scanning Long Press:Enter Remote Control Interface