



# TPC

## Replacement Parts List



Index No.	Part Number	Description
1	L700-0001	Housing cover for transformer
2	P56004-700	110V Power Transformer
3	L700-0002	Transformer Housing
4	L700-0003	Horizontal Arm
5	L700-0004	Counter Sunk Pin
6	L700-0005	Swivel Arm Knuckle
7	L700-0006	Bridge ARM Connector
8	L700-0007	Set Screw (4 X 5 mm)
9	L700-0008	Yoke Joint Cover
10	L700-0009	Balance Arm Assembly
11	L700-0010	Set Screw (5 X 8 mm)
12	L700-0011	Wire Harness Balance Arm
13	L700-0012	Spring For Balance Arm

# TPC



# TPC

Index No.	Part Number	Description
1	L700-703	Light Shield
2	L700-702	Light Handle
3	L700-701	Light Head Yoke
4	L700-704	Frame
5	L700-705	Yoke Cover
6	L700-706	Frame Back Cover
7	L700-707	Yoke Wire Harness
8	L700-708	Power Switch
9	L700-709	Joint Nut
10	L700-710	Flat Washer
11	L700-711	Wave Washer
12	L700-712	Nut For Joint Bolt
13	L700-713	Yoke Side Cover (CAP)
14	L700-714	Spring Holder For Light Bulb
15	L700-P107	Halogen Light Bulb
16	L700-716	Ground Wire
17	L700-717	Lamp Socket
18	L700-718	Socket Housing
19	L700-719	Rubber Washer
20	L700-P705	Reflector
21	L700-721	Nut
22	L700-722	Deflector
23	L700-723	Wire Switch
24	L700-724	Screw
28	L700-723	Rubber She eel For Yoke
29	L700-729	Socket Holder