



# TPC

## Replacement Parts List




Item	Part	Number Description
1	56051	Yoke
2	56052	Support Yoke
3	56053	Handle Frame
4	56054	Socket Bushing
5	56055	Reflector
6	56056	Shield
7	56057	Rear Cover
8	56058	Tension Covers
9	56059	Light Bulb Tension Spring
10	P709	Light Bulb
11	56061	Light bulb Socket
12	56062	Reflector Holder Bushing
13	56063	Washer
14	56064	Rubber Grommet
15	56065	Large Washer
16	56066	Round Nut
17	56067	Deflector
18	56068	Cover, Switch
19	56069	Power Switch
20	56070	Sensor PCB
21	56071	Screw
22	56072	Light Head Wire Harness

# TPC



Item	Part Number	Description
1	56001	Light Bushing
2	56002	Rigid Arm
3	56003	Cover, Transformer
4	56004	Transformer
5	56005	PCB Board
6	56006	Fuse
7	56007	Switch
8	56008	Screw X 4
9	56009	Transformer Cover Screw
10	56010	Set Screw



Item	Part Number	Description
1	56021	Cover Flex Arm
2	56022	Knuckle, Rigid Arm
3	56023	Knuckle, Head
4	56024	Spring Rod
5	56025	Extension Support
6	56026	Bushing
7	56027	Bezel
8	56028	Bushing
9	56029	Block Spring Tensioner
10	56030	Friction Block
11	56031	Tensioner Rod
12	56032	End Cap Cover
13	56033	Retainer Ring
14	56034	Long Pin
15	56035	Retainer Ring Small
16	56036	Short Pin
17	56037	Retainer Ring Small
18	56038	Spring
19	56039	Screw Cap Head Socket
20	56040	Screw, End Cap
21	56041	Friction Screw