

SuperbVac

Wet suction unit



WV6 / WV10



WV3 / WV5



WV3R / WV5R

Features:

- High Power Ring Blower provide Strong and Stable Suction to ensure Continuous Running.
- Copper Impellers Ensure Stable Performance and Longer Running.
- Unique Water Recycle Device to save 75%-80% water.
- Built-in water and air separator can reduce noise level to ensure quiet and silent operation.

SPECIFICATIONS

Model	WV3/WV3R (Single)	WV5/WV5R (Single)	WV6 (Twin)	WV10 (Twin)
Model Water consumption w/o Water Recycling:	4L/min	4L/min	6L/min	6L/min
Water consumption with TPC Water Recycling:	0.8L/min	0.8L/min	1.0L/min	1.0L/min
Suction Rate	900L/min	1500L/min	1800L/min	3000L/min
Vacuum Level (IN Hg)	8-15Hg adjustable, normally is set at 10Hg at the factory			
Noise Level dB(A)	58	58	65	65
Horsepower	1.2HP	2HP	2.4HP	4HP
Voltage	220v	220v	220v	220v
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Max. Users:	3 sets	5 sets	6 sets	10 sets
Current (A)	8A	13A	8A×2	13A×2
Circuit Breaker (A)	20A	20A	20A×2	20A×2
Inlet Connection Size	3/4"	3/4"	1"	1"
Drain Connection Size	32mm	32mm	32mm	32mm
City Water Connection Size	8mm	8mm	8mm	8mm
Air Exhaust Connection Size	60mm	60mm	60mm	60mm
Inlet Water Pressure (PSI)	Above 10PSI or 0.7bar			
Ambient Temperature	2°C —40°C			
Product Size: (L×W×H) mm	450*420*520 580*515*590	450*420*520 580*515*590	720×430×460	720×430×460
Product Weight (kg)	26/31	30/37	60	64

Place Selection Valve

Features:

- Place selection valve with patent technology, easy to install and durable.
- Suitable for any brand of dental chair installation.
- No need to install power supply or solenoid valve.
- Compact size, can be installed in the side box or any space on the ground box.
- Various Suction Tubing Adapters are available.
- Reliable quality with 5 years warranty.

